

Begin your response to **QUESTION 3** on this page.

3. A student is given the task of determining the molar concentration of a  $\text{CuSO}_4$  solution using two different procedures, precipitation and spectrophotometry.

For the precipitation experiment, the student adds 20.0 mL of 0.200  $M$   $\text{Ba}(\text{NO}_3)_2$  to 50.0 mL of the  $\text{CuSO}_4(aq)$ . The reaction goes to completion, and a white precipitate forms. The student filters the precipitate and dries it overnight. The data are given in the following table.

Mass of dry filter paper	0.764 g
Volume of $\text{CuSO}_4(aq)$	50.0 mL
Volume of 0.200 $M$ $\text{Ba}(\text{NO}_3)_2$	20.0 mL
Mass of filter paper and dried precipitate	1.136 g

(a) Write a balanced net ionic equation for the precipitation reaction.

(b) Calculate the number of moles of precipitate formed.

(c) Calculate the molarity of the original  $\text{CuSO}_4$  solution.

**GO ON TO THE NEXT PAGE.**

Use a pencil or pen with black or dark blue ink only. Do NOT write your name. Do NOT write outside the box.

© 2021 College Board.

Visit College Board on the web: collegeboard.org.

Continue your response to **QUESTION 3** on this page.

For the spectrophotometry experiment, the student first makes a standard curve. The student uses a  $0.1000\text{ M}$  solution of  $\text{CuSO}_4(aq)$  to make three more solutions of known concentration ( $0.0500\text{ M}$ ,  $0.0300\text{ M}$ , and  $0.0100\text{ M}$ ) in  $50.00\text{ mL}$  volumetric flasks.

(d) Calculate the volume of  $0.1000\text{ M}$   $\text{CuSO}_4(aq)$  needed to make  $50.00\text{ mL}$  of  $0.0500\text{ M}$   $\text{CuSO}_4(aq)$ .


(e) Briefly describe the procedure the student should follow to make  $50.00\text{ mL}$  of  $0.0500\text{ M}$   $\text{CuSO}_4(aq)$  using  $0.1000\text{ M}$   $\text{CuSO}_4(aq)$ , a  $50.00\text{ mL}$  volumetric flask, and other standard laboratory equipment. Assume that all appropriate safety precautions will be taken.

**GO ON TO THE NEXT PAGE.**

Use a pencil or pen with black or dark blue ink only. Do NOT write your name. Do NOT write outside the box.

Continue your response to **QUESTION 3** on this page.

The standard curve is given below.



(f) The absorbance of the  $\text{CuSO}_4$  solution of unknown concentration is 0.219. Determine the molarity of the solution.

(g) A second student performs the same experiment. There are a few drops of water in the cuvette before the second student adds the  $\text{CuSO}_4(aq)$  solution of unknown concentration. Will this result in a  $\text{CuSO}_4(aq)$  concentration for the unknown that is greater than, less than, or equal to the concentration determined in part (f)? Justify your answer.

**GO ON TO THE NEXT PAGE.**


Use a pencil or pen with black or dark blue ink only. Do NOT write your name. Do NOT write outside the box.

© 2021 College Board.

Visit College Board on the web: collegeboard.org.

**Question 3: Long Answer****10 points**

- (a) For the correct balanced equation (state symbols not required): **1 point**



- (b) For the correct calculated value of the mass of precipitate (may be implicit): **1 point**

$$1.136 \text{ g} - 0.764 \text{ g} = 0.372 \text{ g BaSO}_4$$

For the correct calculated value of the number of moles, consistent with mass of precipitate: **1 point**

$$0.372 \text{ g} \times \frac{1 \text{ mol}}{233.39 \text{ g}} = 0.00159 \text{ mol}$$

**Total for part (b) 2 points**

- (c) For the correct calculated value, consistent with part (b): **1 point**

$$0.00159 \text{ mol BaSO}_4 \times \frac{1 \text{ mol CuSO}_4}{1 \text{ mol BaSO}_4} = 0.00159 \text{ mol CuSO}_4$$

$$\frac{0.00159 \text{ mol CuSO}_4}{0.0500 \text{ L}} = 0.0318 M \text{ CuSO}_4 \quad (0.0319 M \text{ if decimals are carried})$$

- (d) For the correct calculated value: **1 point**

$$M_1V_1 = M_2V_2$$

$$V_1 = \frac{(0.0500 M)(50.00 \text{ mL})}{(0.1000 M)} = 25.0 \text{ mL}$$

- (e) For a correct technique to measure the volume of solution: **1 point**

*First, measure out the correct volume of 0.1000 M CuSO<sub>4</sub> solution with a 25.0 mL volumetric pipet (graduated cylinder or buret is acceptable).*

For a correct technique to dilute the solution to the final volume: **1 point**

*Transfer the 25.0 mL of solution to a 50.00 mL volumetric flask and dilute the solution with water up to the 50.00 mL mark.*

**Total for part (e) 2 points**

- (f) For the correct value (between 0.032 M and 0.038 M): **1 point**

Accept one of the following:

- $y = mx = \frac{0.63}{0.1000}x = 6.3x$

$$x = \frac{y}{6.3} = \frac{0.219 \text{ M}}{6.3} = 0.035 \text{ M}$$

- *Estimated value from the graph within the specified range.*

(g) For the correct answer:	<b>1 point</b>
<i>The concentration will be <u>less than</u> that determined in part (f).</i>	
For a valid justification:	<b>1 point</b>
<i>The additional water will decrease the concentration of CuSO<sub>4</sub> in the cuvette. Therefore, there will be a decrease in absorbance (according to the Beer-Lambert law). This dilution results in a lower estimated concentration of CuSO<sub>4</sub>.</i>	
<b>Total for part (g)</b>	<b>2 points</b>
<b>Total for question 3</b>	<b>10 points</b>